



| O2k and DatLab file: P___ (A / B) | | 2017- | | Operator: | | | | |
|--|------------------------------|----------------------------|-------------------------|----------------|---------------------------|--------------|---|---|
| Sample type: | | Cohort: | | Sample code: | | | | |
| Sample.Subsample number: | | Unit: | | Concentration: | | | | |
| Medium: MiR05Cr | | | | | | | | |
| Event | Mark name O2/Amp Slope | Mark name Amp Signal | Final conc. 2 ml O2k | Stock [mM] | Comment | Tit. [µl] | A | B |
| | | | ~200 µM O2 | | MiR | | | |
| 0AmR | | | 10 µM | 10 | | 2 | | |
| 0HRP | | | 1 U/ml | | 500 U/ml stock | 4 | | |
| 0SOD | 0BG | | 5 U/ml | | 5000 U/ml stock | 2 | | |
| (0.0) | | 0.0 | 0 µM HP | | <i>Mark on Amp signal</i> | | | |
| (0.1) | | 0.1 | 0.1 µM HP | 0.04 | <i>Mark on Amp signal</i> | 5 | | |
| (0.2) | | 0.2 | 0.2 µM HP | 0.04 | <i>Mark on Amp signal</i> | 5 | | |
| 1ce | 1R | | | | | | | |
| | | | 0 | | | | | |
| 1Dig | | | | 8.1 | | | | |
| (1.0) | | 1.0 | 0 µM HP | | <i>Mark on Amp signal</i> | | | |
| (1.1) | | 1.1 | 0.1 µM HP | 0.04 | <i>Mark on Amp signal</i> | 5 | | |
| 1P | | | 5 mM | 2000 | | 5 | | |
| 1M | 1PM | | 2 mM | 400 | | 10 | | |
| 2D | 2D | | 1/2.5 mM | 500 | | 4/10 | | |
| 3U* | 3U | | Δ0.5 µM | 1 | CCCP | Δ1 | | |
| 4G | 4G | | 10 mM | 2000 | | 10 | | |
| (2.0) | | 2.0 | 0 µM HP | | <i>Mark on Amp signal</i> | | | |
| (2.1) | | 2.1 | 0.1 µM HP | 0.04 | <i>Mark on Amp signal</i> | 5 | | |
| 5S | 5S | | 50 mM | 1000 | | 100 | | |
| 6Oct | 6Oct | | 0.5 mM | 100 | | 10 | | |
| 7Rot | 7Rot | | 0.5 µM | 1 | | 1 | | |
| (3.0) | | 3.0 | 0 µM HP | | <i>Mark on Amp signal</i> | | | |
| (3.1) | | 3.1 | 0.1 µM HP | 0.04 | <i>Mark on Amp signal</i> | 5 | | |
| 8Gp | 8Gp | | 10 mM | 1000 | | 20 | | |
| 9Ama | 9Ama | | 2.5 µM | 5 | | 1 | | |
| | | 5HPO | 0 | | <i>Mark on Amp signal</i> | | | |
| 5HP.1 | | 5HP.1 | 0.1 µM | 0.04 | <i>Mark on Amp signal</i> | 5 | | |