Technical Communication

Magnesium Green for fluorometric measurement of ATP production does not interfere with

4 mitochondrial respiration

- 5 Duiza HD Cardoso*, Carolina Doerrier, Erich Gnaiger
- 6 Oroboros Instruments, Innsbruck, Austria 7
- 8 *Corresponding author: luiza.cardoso@oroboros.at

Abstract

1

9 10

20 21

22

23

24

25

26

27

28

29

30

31 32

33

34

35 36

37

38

39 40

- For the advanced study of mitochondrial function, high-resolution respirometry is 11 extended by fluorometric measurement of ATP production using the fluorophore 12 Magnesium Green™ (MgG). A common problem with several fluorescent dyes is the 13 inhibition of mitochondrial respiration. In the present study, a coupling control 14 15 protocol was applied in combination with MgG to measure ATP production simultaneously with respiration for calculation of P»/O2 ratios. MgG at 1.1 µM did 16 not affect respiration through the NADH-linked and succinate-linked pathways. 17 Respiration was not inhibited in any of the coupling control states, hence coupling 18 19 control efficiencies were not affected by MgG.
 - *Keywords* ATP; ATP production; high-resolution respirometry; Magnesium Green; mitochondria; oxidative phosphorylation; fluorometry; FluoRespirometry

1. Introduction

Mitochondrial ATP production can be analyzed with a fluorometric technique using Magnesium Green™ (MgG) as a fluorescent probe, as described by Chinopoulos et al (2009). Application of the Mg²+-sensitive fluorophore as an indicator of ATP production relies on the fact that ADP and ATP have different affinities for Mg²+ (Gnaiger, Wyss 1994; Leyssens et al 1996; Budinger et al 1998). ADP is phosphorylated to ATP in the mitochondrial matrix. In the phosphorylation system ADP/ATP and inorganic phosphate Pi are exchanged stoichiometrically by the adenine nucleotide translocase ANT and the phosphate carrier PiC. Under experimental conditions when ADP decreases while ATP increases in the extramitochondrial milieu, the Mg²+ concentration declines due to the higher affinity for Mg²+ of ATP than ADP (Figure 1). Therefore, the fluorometric assay with the membrane-impermeant MgG provides a quantitative approach to analyze mitochondrial ATP production. This method was developed further to measure concomitantly mitochondrial ATP production and O₂ consumption in the Oroboros O2k-FluoRespirometer which is an experimental system complete for high-resolution respirometry including fluorometry (Chinopoulos et al 2014).

www.mitofit.org 1 of 12

Fluorescent dyes are widely used to assess various parameters relevant in mitochondrial physiology. Safranin, rhodamine and its derivatives, such as TMRM, are frequently employed as reporters of the mitochondrial membrane potential $\Delta\Psi_{p+}$. However, all $\Delta\Psi_{p+}$ dyes have been shown to affect mitochondrial respiration (Scaduto, Grotyohann 1999). Like TPP+, safranin mainly affects the NADH (N)-linked pathway, the phosphorylation system, and to a smaller extent the succinate (S)-linked pathway (Krumschnabel et al 2014). The effect of $\Delta\Psi_{p+}$ fluorescent probes can be explained since they accumulate in the mitochondrial matrix and thus possibly affect mitochondrial function.




Figure 1. Concept of the MgG assay according to Chinopoulos et al (2014). MgG fluoresces when bound to Mg²⁺. ADP and ATP compete for Mg²⁺ binding with different affinities; ATP has a higher affinity for Mg²⁺ compared to ADP. (a) Initial experimental conditions, when ADP is added, binding some Mg²⁺, and high MgG fluorescence drops slightly. (b) As the experiment proceeds, ADP is phosphorylated to ATP, which is exchanged for ADP by the adenine nucleotide translocase ANT. With increase in ATP in the extramitochondrial medium, more Mg²⁺ is bound to ATP, and MgG fluorescence decreases.

It is important to note that another dye frequently used in mitochondrial physiology studies, Amplex UltraRed, employed to analyze $\rm H_2O_2$ production, was shown to affect mitochondrial respiration even though the mitochondrial membranes are not permeable to this fluorophore (Makrecka-Kuka et al 2015). Therefore, it is important to analyze whether MgG affects mitochondrial respiration, despite the fact that mitochondrial membranes are not permeable to this fluorophore.

In the present technical communication, we report the effect of MgG on mitochondrial respiration, which is the gold standard to evaluate mitochondrial function. This provides an important contribution towards further development of this method to analyze $P \gg /O_2$ ratios in different mitochondrial preparations. The use of a coupling control protocol assessing O_2 consumption and MgG fluorescence allows for the evaluation of mitochondrial respiration and ATP production using NADH- and succinate-linked substrates in LEAK, OXPHOS- and ET-state, making it possible to obtain flux control efficiencies.

2. Materials and methods

2.1. Reagents

75 76 77

78

79

80

81

82

83

84 85

86

87

88

89

90 91 92

93

94

95

96 97

98 99

100

101 102

103 104

105

106

107

108

109

110

111

112

113

114

115

116

117

118

119

120 121 Magnesium Green was purchased from Invitrogen/Thermo Fisher Scientific (cat. Nº M3733). Antimycin A (cat. Nº A8674), ATP (cat. Nº A2383), CCCP (cat. Nº C2759), malate (cat. Nº M1000), MgCl₂ 1 M (cat. Nº M1028), oligomycin (cat. Nº O4876), pyruvate (cat. Nº P2256), rotenone (cat. Nº R8875), SF 6847 (cat. Nº T182), and succinate (cat. Nº S2378) were obtained from Sigma Aldrich. ADP was acquired from Millipore (cat. Nº 117105), and carboxyatractyloside from Calbiochem (cat. Nº 216201).

ADP and ATP were diluted in deionized H_2O without addition of Mg^{2+} salts, pH was adjusted to 6.9 with KOH. Magnesium Green, malate, succinate, carboxyatractyloside and $MgCl_2$ were diluted in deionized H_2O whereas antimycin A, CCCP, oligomycin, rotenone and SF 6847 were diluted in ethanol p.a. All solutions were aliquoted and stored at -20 °C, except pyruvate, which was diluted in deionized H_2O fresh on the day of each experiment.

2.2. Animals

Wild-type C57BL/6N adult mice (*N*=3 per experimental group) were housed in the animal facility of the Medical University of Innsbruck (maximum 5 mice per cage) and, maintained at 22 °C with a controlled 12 h light/dark cycle. Mice were fed *ad libitum* with free access to water. All procedures were conducted according to the Austrian Animal Experimentation Act in compliance with the European convention for the protection of vertebrate animals used for experimental and other scientific (Tierversuchsgesetz 2012; Directive 2010/63/EU: BMWFM-66.011/0128-WF/V/3b/2016).

2.3. Cardiac mitochondrial isolation and protein concentration determination

Following cervical dislocation, the hearts were immediately excised and transferred into ice-cold biopsy preservation solution (BIOPS: 10 mM Ca²⁺-EGTA - 0.1 μM free Ca²⁺, 20 mM imidazole, 20 mM taurine, 50 mM K+-MES, 0.5 mM dithiothreitol, 6.56 mM MgCl₂, 5.77 mM ATP, 15 mM phosphocreatine, pH 7.1 adjusted with KOH) for short period of time (1-2 h; Fontana-Ayoub et al 2016). All procedures were performed on ice (Gnaiger et al 2000a). Mouse heart mitochondria were isolated following the protocol described by Fontana-Ayoub and Krumschnabel (2015). The heart (~ 80–120 mg) was washed to remove blood clots and minced with 1 mL of BIOPS. The tissue was homogenized with 2 mL isolation buffer (IB1: 0.5 M mannitol; 0.5 M sucrose; 0.1 M EGTA; pH 7.4 adjusted with Tris; 2.5 mg/mL BSA and 0.5 mg/mL subtilisin, the latter two added freshly on the day of use) on a 10 mL glass-Teflon Potter Elvehjem homogenizer, 6-8 × with about 1000 rpm mechanical rotation. 3 mL of IB1 was added to the homogenate which was centrifuged at 800 g for 10 min at 4 °C. The supernatant was centrifuged again, at 10 000 g for 10 min at 4 °C. The pellet was resuspended carefully using a 1 mL pipette in 0.5 mL IB2 (IB1 without subtilisin). After addition of 2 mL IB2, the homogenate was centrifuged again at 10 000 g for 10 min at 4 °C. The pellet was resuspended in 200 μL of IB3 (IB1 without BSA and subtilisin) and kept on ice until use on the same day within 2 h.

www.mitofit.org 3 of 12



Protein concentration was used for calculation of mass specific O_2 flux, determined using the kit DC Protein Assay (Bio-Rad, Hercules, CA, US). Absorbance was measured at 620 nm with a Tecan Infinite TM F200 spectrophotometer (Tecan, Männedorf, Switzerland), using BSA at different concentrations as standards (Lowry et al 1951).

2.4. High-resolution respirometry

Oxygen consumption and ATP production measurements were performed simultaneously at 37 °C in the O2k-FluoRespirometer (O2k, Oroboros Instruments, Innsbruck, Austria). The O2k includes two Duran® glass chambers with stirring (750 rpm) and controlled temperature for closed-chamber respirometry using polarographic oxygen sensors (POS). Smart Fluo-Sensors Blue were used, with excitation LED 465 nm and filters for the LED and photodiode selected for Magnesium GreenTM). Specific amperometric emission and detection settings — fluorescence light intensity of 500 and gain 100 — were applied with the software DatLab 7.4 (Oroboros Instruments, Innsbruck, Austria) with continuous data recording set at 2 s time intervals. Standardized calibrations and instrumental O_2 background tests were performed (Doerrier et al 2018). The time-derivative of the O_2 concentration is calculated real-time by DatLab, providing traces of O_2 flux corrected for the O_2 instrumental background (Gnaiger 2001).

Experiments were run with cardiac isolated mitochondria at protein concentrations in the range of 0.026-0.049~mg/mL in modified mitochondrial respiration medium MiR05-MgG (MgCl2 1 mM instead of 3 mM in MiR05, EGTA 0,5 mM, KH₂PO₄ 10 mM, Hepes 20 mM, lactobionic acid 60 mM, D-sucrose 110 mM, taurine 20 mM, BSA 1 g/L, pH adjusted with KOH to 7.1). This modification of MiR05 (Gnaiger et al 2000a) was optimized for measurement of ATP production with MgG.

2.5. ATP production measurement with MgG

MgG (Magnesium GreenTM, pentapotassium salt, cell impermeant) does not permeate biological membranes. Therefore, the plasma membrane barrier function must be removed, as achieved in mitochondrial preparations – isolated mitochondria, tissue homogenates, permeabilized tissues and cells. MgG remains outside of the mitochondrial matrix and fluoresces when bound to Mg^{2+} . In the phosphorylation reaction

$$ADP + P_i \leftrightarrows ATP$$

reactants and MgG bind Mg²⁺ according to their apparent dissociation constants. When ADP is added to the experimental chamber, there is a fast drop of the fluorescence signal. If mitochondria and fuel substrates are present, ATP is generated and exchanged with ADP by the ANT. ATP has a higher affinity to Mg²⁺ compared to ADP. As ATP concentration increases in the medium, the free Mg²⁺ concentration declines, less MgG is bound to Mg²⁺, and the fluorescence decreases. The ATP concentration in the medium is calculated according to Chinopoulos et al (2009; 2014), taking in account that: (1) the initial concentration of ATP is zero, (2) the initial concentration of ADP is known, (3) the concentration of Mg²⁺ is measured, and (4) apparent K_d values for ADP and ATP with Mg²⁺ are obtained experimentally.



The free Mg^{2+} concentration was calibrated in MiR05-MgG containing the mitochondrial sample, fuel substrates, carboxyatractyloside, and oligomycin. MgCl₂ was titrated in 10 steps of 0.1 mM to obtain a non-linear fit for calibration of the amperomeric signal. After calibration, the K_d of ADP and ATP for Mg^{2+} was determined for each experimental condition by performing multiple titrations with ADP or ATP.

2.6. Substrate-uncoupler-inhibitor-titration (SUIT) protocols

Coupling control protocols (SUIT-006) assess different coupling control states -LEAK, OXPHOS and ET - at a constant electron-transfer-pathway state (Gnaiger et al 2020). The effect of MgG on mitochondrial respiration was evaluated in its absence or presence (1.1 µM), which was added to the experimental chambers prior to sample addition. Since this fluorescent dye is diluted in water, and only a 2 µL volume was added into the 2 mL chamber, no solvent addition was performed in the control group without MgG. After addition of isolated mitochondria into the O2k chambers, residual oxygen consumption Rox was measured in the absence of substrates. Two coupling-control protocols were used to study simultaneously oxygen consumption and ATP production with the following titrations: NADH-pathway with 5 mM pyruvate and 2 mM malate, or Succinate-pathway with 0.5 µM rotenone and 10 mM succinate. First, LEAK respiration was measured in the absence of ADP. Secondly, OXPHOS capacity was measured after addition of 2 mM ADP. Oligomycin (7.5–10.0 nM) or carboxyatractyloside (0.3 – 0.4 mM) were added to induce again a LEAK state. This was followed by stepwise titration of the uncouplers CCCP (0.5 µM steps) or SF 6847 (25-50 nM steps) up to the optimum concentration, when the maximum O₂ flux was achieved as a measure of ET capacity. Finally, residual oxygen consumption was measured after the addition of the CIII inhibitor antimycin A (2.5 μ M).

2.7. Data analysis

The assays were repeated 3 times with independent mitochondrial preparations, with or without MgG, for each condition tested. Data analysis for O_2 consumption, calculations of K_d values and ATP production following Chinopoulos et al 2014, were performed using the templates provided with the software DatLab 7.4.

3. Results and discussion

Figures 2a and 3a show superimposed traces of O_2 concentration and O_2 flux per mass. Coupling control of mitochondrial respiration was measured in two different electron-transfer-pathway control states. In the N-protocol, the NADH-linked pathway through Complex I (CI) was evaluated in the presence of pyruvate and malate which stimulate dehydrogenases of the TCA cycle, leading to reduction of NAD+ to NADH. NADH is the substrate of CI, with further electron flow into the Q-junction, CIII and CIV (Figure 2). In the S-protocol, CI was inhibited by rotenone to prevent reverse electron transfer and accumulation of oxaloacetate, which is an inhibitor of succinate dehydrogenase (Makrecka-Kuka et al 2015; Gnaiger 2020), and respiration was measured supported by succinate as the substrate of CII (Figure 3).

www.mitofit.org 5 of 12

In both protocols, LEAK respiration was measured (1) L(n), in the absence of adenylates and (2) L(Omy) or L(Cat), in the presence of phosphorylation system inhibitors. Respiration in these two LEAK states was similar, but slightly lower in L(Omy) with the N-protocol (Figure 2a, Table 1). L(n) stabilized quickly, whereas for L(Omy) it took a long time to fully inhibit respiration by the low concentration of 7.5–10.0 nM oligomycin. In the S-protocol with sequential addition of rotenone followed by succinate, L(n) increased for a few minutes until stabilization (Figure 3a). Inhibition by carboxyatractyloside (0.3–0.4 μ M) was immediate, and L(Cat) tended to be slightly lower than L(n) (Table 1).




Figure 2. Simultaneous measurement of respiration and ATP production by high-resolution FluoRespirometry in mitochondria isolated from mouse heart. Representative traces for coupling control protocol SUIT-006 with NADH-linked substrates (N-protocol), following additions (respiratory states): isolated mitochondria imt (ROX), pyruvate & malate (LEAK), ADP (OXPHOS), oligomycin (LEAK), uncoupler U (ET), and antimycin A (ROX). Experiment 2019-02-07 P5 04: (a) O_2 concentration (dark and lighter blue traces) and O_2 flux per mass (dark and lighter red), 1.1 μ M MgG versus control; (b) MgG fluorescence signal; (c) ATP concentration calculated from MgG signal calibrated as Mg²⁺ concentration.




Figure 3. Simultaneous measurement of respiration and ATP production by high-resolution FluoRespirometry in mitochondria isolated from mouse heart. Representative traces for coupling control protocol SUIT-006 with succinate as substrate (S-protocol), following additions (respiratory states): isolated mitochondria imt (ROX), succinate & rotenone (LEAK), ADP (OXPHOS), carboxyatractyloside (LEAK), uncoupler U (ET), and antimycin A (ROX). Experiment 2019-03-18 P5-03: (a) O_2 concentration (dark and lighter blue traces) and O_2 flux per mass (dark and lighter red), 1.1 μ M MgG versus control; (b) MgG fluorescence signal; (c) ATP concentration calculated from MgG signal calibrated as Mg²⁺ concentration.

OXPHOS capacity P was measured in the presence of a kinetically saturating concentration of ADP. The optimum uncoupler concentrations to measure maximum ET capacity E were 6.0–7.0 μ M CCCP in the N-protocol, and 0.150–0.175 μ M SF 6847 in the S-protocol. In the N-protocol, P was stable over time and identical to E. However, in the S-protocol, P showed a slight decrease over time. Extrapolating this trend of declining O_2 flux to the point where ET capacity was measured explains why E appears to be lower than P (dashed trendline, Figure 3a). In both protocols, therefore, E = P, indicating that OXPHOS capacity was not limited by the phosphorylation system. This agrees with results

www.mitofit.org 7 of 12



 for mouse heart mitochondria on coupling control even in the combined NS-pathway (Lemieux et al 2017). Parallel measurements were performed in the presence and absence of 1.1 μ M MgG with the N- and S-protocol. This low concentration of MgG used is sufficiently high for calculating ATP production (Figures 2b and c and Figures 3b and c).


The MgG assay to measure ATP production can be used concomitantly with high-resolution respirometry, providing information real-time. Other methods are available to detect ATP production real-time. Spectrophotometric detection of NADPH can be used in conjunction with the coupled enzyme system hexokinase and glucose-6-phosphate dehydrogenase (Horgan, 1978). This assay has been adapted for simultaneous detection of O₂ consumption and NADPH (Lark et al 2016). The luciferin/luciferase assay can be used for continuous measurement of ATP production (Manfredi et al 2002). It is important to note that luciferase consumes O₂, and instruments typically used for luminometry do not allow monitoring of O₂ concentration in parallel.

Another method for continuous measurement of the $P*/O_2$ ratio is the steady-state ADP injection-respirometry (Gnaiger et al 2000b; 2001). The phosphorylation rate is set by continuous injection of ADP as the rate-limiting step while measuring O_2 consumption stimulated to a constant sub-maximal level. Chance and Williams (1955) originally described a polarographic ADP pulse-titration method to determine the $P*/O_2$ ratio, titrating a known concentration of ADP, which leads to a peak of O_2 consumption stimulated by the complete phosphorylation of ADP to ATP. The ADP pulse-titration method has been extended and critically discussed by Gnaiger (2001).

End-point assays are available to detect ATP levels, providing discontinuous measurement of ATP production. These include chromatography (high performance liquid chromatography, HPLC; thin layer chromatography, TLC); nuclear magnetic resonance detection of 2-deoxyglucose and its phosphorylated form, and radioactivity measurements using ³²P (Menegollo et al 2019; Morciano et al 2017; Fink et al 2017; Sausen et al 2019).

The fluorometric MgG assay applied simultaneously with O_2 consumption by HRR has been used extensively (Iftikar, Hickey 2013; Goo et al 2013; Chinopoulos et al 2014; Pham et al 2014; Power et al 2014; Salin et al 2016; Napa et al 2017; Masson et al 2017; Salin et al 2018; Devaux et al 2019; Salin et al 2019). Understanding whether MgG may affect respiration is crucial for such studies, particularly for $P_{\rm w}/O_2$ ratios obtained in different electron-transfer-pathway states.

It is well established that different dyes commonly applied to measure mitochondrial membrane potential inhibit OXPHOS capacity, *e.g.*, safranin, rhodamine 123 and its derivatives TMRM and TMRE (Krumschnabel et al 2014; Scaduto, Grotyohann 1999). Surprisingly, Amplex UltraRed used to detect $\rm H_2O_2$ flux impairs respiration despite not accumulating in the mitochondria (Makrecka-Kuka et al 2015). Therefore, we studied the effect of MgG on respiration. MgG at 1.1 μM did not affect NADH-linked nor succinate-linked respiration in any coupling control state (LEAK, OXPHOS and ET) measured in mitochondria isolated from mouse hearts (Figure 4). In addition, residual oxygen consumption was not affected by MgG.






Figure 4. O² **consumption in the absence and presence of MgG by mitochondria isolated from mouse heart.** The respiratory rates indicated in the abscissa were measured by HRR with two coupling control protocols SUIT-006, with the following respiratory states: ROX, LEAK (in the absence of adenylates), OXPHOS, LEAK (in the presence of inhibitors), ET, and ROX. Sequential titrations are described for **(a)** N-protocol (experiments 2019-02-05 P3-04, 2019-02-06 P3-03 and 2019-02-07 P5-04) and **(b)** S-protocol (experiments 2019-03-13 P6-03, 2019-03-14 P3-03 and 2019-03-18 P5-03). For both graphs the three symbol shapes show independent mitochondrial preparations, whereas open and closed symbols compare results in controls and in the presence of MgG from the same preparation; bars represent the average.

www.mitofit.org 9 of 12



Table 1. Coupling control efficiency (P-L)/P and $P \gg /O_2$ ratio in absence or presence of MgG. Average \pm SD, N=3. OXPHOS capacity P and LEAK respiration L corrected for residual oxygen consumption Rox. L(n)/L(inh) ratios: L in the absence of adenylates (n) over L with an inhibitor (inh) of the phosphorylation system, oligomycin Omy or carboxyatractyloside Cat for the N- or S-pathway, respectively. L(inh) is used in (P-L)/P.

Protocol	(P-L)/P	L(n)/L(inh)	P»/O ₂	P»/0
N-pathway - MgG	0.90 ± 0.01	1.13 ± 0.05	-	-
N-pathway + MgG	0.88 ± 0.02	1.12 ± 0.04	2.33 ± 1.07	1.16 ± 0.53
S-pathway - MgG	0.62 ± 0.05	1.27 ± 0.16	-	-
S-pathway + MgG	0.61 ± 0.05	1.12 ± 0.26	2.78 ± 0.74	1.39 ± 0.37

The NADH-pathway has three coupling sites, CI, CIII and CIV, whereas the succinate-pathway has only the latter two, resulting in a lower P_{P}/O_{2} ratio. When dividing ATP flux, calculated from the increase in ATP concentration per time, by the simultaneously measured O_{2} flux, then P_{P}/O_{2} flux ratios $(J_{P}/J_{O_{2}})$ are obtained. The P_{P}/O_{2} is twice the classical P_{P}/O_{2} (Table 1). P_{P}/O_{2} obtained for S-pathway was close to the theoretically expected value (Gnaiger et al 2020). The result obtained for N-pathway was lower than expected. A limitation of the present study is the low number of replicates (N = 3), with a high variability of P_{P}/O_{2} ratios. Further experiments are in preparation.

Coupling control efficiencies are closely related to $P*/O_2$ ratios. The coupling control efficiency is defined as (E-L)/E, ranging from 0, at zero coupling, to 1 in a fully coupled system. In the present case of P=E, the coupling control efficiency is expressed as the P-L control efficiency, (P-L)/P (Gnaiger 2020). As expected, a higher P-L control efficiency of 0.89 ± 0.02 was found for the N-pathway than 0.62 ± 0.05 for the S-pathway (pooled data with and without MgG, average \pm standard deviation, N=6; Table 1). These correspond to a RCR = P/L of 9.6 ± 1.8 for the N-pathway and 2.6 ± 0.3 for the S-pathway.

In summary, MgG did not affect respiration in any of the coupling control states. These results demonstrate that measurement of O_2 consumption is reliable concomitant with the MgG assay in SUIT protocols with different pathway states and coupling states.

Acknowledgements

We thank Marco Di Marcello and Manuela Passrugger for expert technical support on media and chemicals preparation, equipment maintenance and mitochondria isolation. This work was partially funded by the European Union's Horizon 2020 research and innovation programme under grant agreement No. 859770, NextGen-O2k project. Contribution to COST Action CA15203 MitoEAGLE.

Author contributions

LHDC, CD and EG designed the work; LHDC collected and analyzed data and drafted the article; CD and EG critically revised the article, all authors approved the final version of the manuscript.

Conflicts of interest

EG is founder and CEO of Oroboros Instruments, Innsbruck, Austria.

Data availability

Original files are available Open Access at Zenodo repository: 10.5281/zenodo.4032674.



Abbreviations

Amp amperometric; ANT adenosine nucleotide translocase; BSA bovine serum albumin; CI to CIV Complex I to IV; CCCP carbonyl cyanide m-chlorophenyl hydrazone; $\Delta\Psi_p+$ mt-membrane potential; EGTA ethylene glycol tetraacetic acid; E ET capacity; ETS electron transfer system; F_0F_1 ATP synthase; Hepes N-(2-hydroxyethyl)piperazine-N'-(2-ethanesulfonic acid); HRR high-resolution respirometry; imt isolated mitochondria; J_{02} O_2 flux; K_d dissociation constant; L LEAK respiration; LED light-emitting diode; MES 2-(N-morpholino)ethanesulfonic acid hydrate; MgG Magnesium Green; P OXPHOS capacity; P-N0 ADP phosphorylated per atom oxygen consumed; P-N1 inorganic phosphate; RCR respiratory acceptor control ratio; R1 R2 R3 residual oxygen consumption; SUIT substrate-uncoupler-inhibitor-titration; TCA tricarboxylic acid; TMRM tetramethylrhodamine methyl ester; TMRE tetramethylrhodamine ethyl ester; TPP+ tetraphenylphosphonium; Tris 2-amino-2-(hydroxymethyl)-1,3-propanediol; U uncoupler.

References

- Budinger GRS, Duranteau J, Chandel NS, Schumacker PT (1998) Hibernation during hypoxia in cardiomyocytes. Role of mitochondria as the O_2 sensor. J Biol Chem 273:3320-6.
- Chance B, Williams GR (1955) Respiratory enzymes in oxidative phosphorylation. I. Kinetics of oxygen utilization. J Biol Chem 217:383-93.
- Chinopoulos C, Vajda S, Csanady L, Mandi M, Mathe K, Adam-Vizi V (2009) A Novel Kinetic Assay of Mitochondrial ATP-ADP Exchange Rate Mediated by the ANT. Biophys J 96:2490-504.
- Chinopoulos C, Kiss G, Kawamata H, Starkov AA (2014) Measurement of ADP-ATP exchange in relation to mitochondrial transmembrane potential and oxygen consumption. Methods Enzymol 542:333-48.
- Devaux JBL, Hedges CP, Birch N, Herbert N, Renshaw GMC, Hickey AJR (2019) Acidosis maintains the function of brain mitochondria in hypoxia-tolerant triplefin fish: a strategy to survive acute hypoxic exposure? Front Physiol 9:1941.
- Doerrier C, Garcia-Souza LF, Krumschnabel G, Wohlfarter Y, Mészáros AT, Gnaiger E (2018) High-Resolution FluoRespirometry and OXPHOS protocols for human cells, permeabilized fibers from small biopsies of muscle, and isolated mitochondria. Methods Mol Biol 1782:31-70.
- Fink BD, Bai F, Yu L, Sivitz WI (2017) Regulation of ATP production: dependence on calcium concentration and respiratory state. Am J Physiol Cell Physiol 313:C146-53.
- Fontana M, Krumschnabel G (2015) Isolation of mouse heart mitochondria. Mitochondr Physiol Network 20.06(01):1-2.
- Fontana-Ayoub M, Fasching M, Gnaiger E (2016) Selected media and chemicals for respirometry with mitochondrial preparations. Mitochondr Physiol Network 03.02(18):1-10.
- Gnaiger E, Wyss M (1994) Chemical forces in the cell: Calculation for the ATP system. In: What is Controlling Life? (Gnaiger E, Gellerich FN, Wyss M, eds) Modern Trends in BioThermoKinetics 3. Innsbruck Univ Press:207-12.
- Gnaiger E (2001) Bioenergetics at low oxygen: dependence of respiration and phosphorylation on oxygen and adenosine diphosphate supply. Respir Physiol 128:277-97.
- Gnaiger E (2020) Mitochondrial pathways and respiratory control. An introduction to OXPHOS analysis. 5th ed. Bioenerg Commun 2020.2: 112 pp. doi:10.26124/bec:2020-0002.
- Gnaiger E et al MitoEAGLE Task Group (2020) Mitochondrial physiology. Bioenerg Commun 2020.1. doi:10.26124/bec:2020-0001.v1.
 - Gnaiger E, Kuznetsov AV, Schneeberger S, Seiler R, Brandacher G, Steurer W, Margreiter R (2000a) Mitochondria in the cold. In: Life in the Cold (Heldmaier G, Klingenspor M, eds) Springer, Heidelberg, Berlin. New York: 431-42.
- Gnaiger E, Méndez G, Hand SC (2000b) High phosphorylation efficiency and depression of uncoupled respiration in mitochondria under hypoxia. Proc Natl Acad Sci U S A 97:11080-5.
- Goo S, Pham T, Han JC, Nielsen P, Taberner A, Hickey A, Loiselle D (2013) Multiscale measurement of cardiac energetics. Clin Exp Pharmacol Physiol 40:671-81.

www.mitofit.org 11 of 12



- Horgan DJ (1978) A spectrophotometric assay of ATP synthesized by sarcoplasmic reticulum. Aust J Biol Sci 31:21-4.
- Iftikar FI, Hickey AJ (2013) Do mitochondria limit hot fish hearts? Understanding the role of mitochondrial function with heat stress in Notolabrus celidotus. PLoS One 8:e64120.
 - Krumschnabel G, Eigentler A, Fasching M, Gnaiger E (2014) Use of safranin for the assessment of mitochondrial membrane potential by high-resolution respirometry and fluorometry. Methods Enzymol 542:163-81.
 - Lark DS, Torres MJ, Lin CT, Ryan TE, Anderson EJ, Neufer PD (2016) Direct real-time quantification of mitochondrial oxidative phosphorylation efficiency in permeabilized skeletal muscle myofibers. Am J Physiol Cell Physiol 311:C239-45.
 - Lemieux H, Blier PU, Gnaiger E (2017) Remodeling pathway control of mitochondrial respiratory capacity by temperature in mouse heart: electron flow through the Q-junction in permeabilized fibers. Sci Rep 7:2840, DOI:10.1038/s41598-017-02789-8.
 - Leyssens A, Nowicky AV, Patterson L, Crompton M, Duchen MR (1996) The relationship between mitochondrial state, ATP hydrolysis, [Mg²⁺]_i and [Ca²⁺]_i studied in isolated rat cardiomyocytes. J Physiol 496:111-28.
 - Lowry OH, Rosebrough NJ, Farr AL, Randall RJ (1951) Protein measurement with the Folin phenol reagent. J Biol Chem 193:265-275.
 - Makrecka-Kuka M, Krumschnabel G, Gnaiger E (2015) High-resolution respirometry for simultaneous measurement of oxygen and hydrogen peroxide fluxes in permeabilized cells, tissue homogenate and isolated mitochondria. Biomolecules 5:1319-38.
 - Manfredi G, Yang L, Gajewski CD, Mattiazzi M (2002) Measurements of ATP in mammalian cells. Methods 26:317-26.
 - Masson SWC, Hedges CP, Devaux JBL, James CS, Hickey AJR (2017) Mitochondrial glycerol 3-phosphate facilitates bumblebee pre-flight thermogenesis. Sci Rep 7:13107.
 - Menegollo M, Tessari I, Bubacco L, Szabadkai G (2019) Determination of ATP, ADP, and AMP Levels by Reversed-Phase High-Performance Liquid Chromatography in Cultured Cells. Methods Mol Biol 1925:223-32.
 - Morciano G, Sarti AC, Marchi S, Missiroli S, Falzoni S, Raffaghello L, Pistoia V, Giorgi C, Di Virgilio F, Pinton P (2017) Use of luciferase probes to measure ATP in living cells and animals. Nat Protoc 12:1542-62.
 - Napa K, Baeder AC, Witt JE, Rayburn ST, Miller MG, Dallon BW, Gibbs JL, Wilcox SH, Winden DR, Smith JH, Reynolds PR, Bikman BT (2017) LPS from *P. gingivalis* negatively alters gingival cell mitochondrial bioenergetics. Int J Dent 2017:2697210.
 - Pham T, Loiselle D, Power A, Hickey AJ (2014) Mitochondrial inefficiencies and anoxic ATP hydrolysis capacities in diabetic rat heart. Am J Physiol 307:C499–507.
 - Power A, Pearson N, Pham T, Cheung C, Phillips A, Hickey A (2014) Uncoupling of oxidative phosphorylation and ATP synthase reversal within the hyperthermic heart. Physiol Rep pii:e12138.
 - Salin K, Villasevil EM, Auer SK, Anderson GJ, Selman C, Metcalfe NB, Chinopoulos C (2016) Simultaneous measurement of mitochondrial respiration and ATP production in tissue homogenates and calculation of effective P/O ratios. Physiol Rep 10.14814/phy2.13007.
 - Salin K, Villasevil EM, Anderson GJ, Selman C, Chinopoulos C, Metcalfe NB (2018) The RCR and ATP/O indices can give contradictory messages about mitochondrial efficiency. Integr Comp Biol 58:486-94.
 - Salin K, Villasevil EM, Anderson GJ, Lamarre SG, Melanson CA, McCarthy I, Selman C, Metcalfe NB (2019) Differences in mitochondrial efficiency explain individual variation in growth performance. Proc Biol Sci 286:20191466.
 - Sausen CW, Rogers CM, Bochman ML (2019) Thin-Layer Chromatography and Real-Time Coupled Assays to Measure ATP Hydrolysis. Methods Mol Biol 1999:245-253.
 - Scaduto RC Jr, Grotyohann LW (1999) Measurement of mitochondrial membrane potential using fluorescent rhodamine derivatives. Biophys J 76:469-77.

Copyright: © 2021 The authors. This is an Open Access preprint (not peer-reviewed) distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original authors and source are credited. © remains with the authors, who have granted MitoFit Preprints an Open Access publication license in perpetuity.

