



O2k and DatLab file: P___ (A / B)		2017-	Operator:					
Sample type:		Cohort:	Sample code:					
Sample.Subsample number:		Unit:	Concentration:					
Medium: MiR05Cr								
Event	Mark name O2/Amp Slope	Mark name Amp Signal	Final conc. 2 ml O2k	Stock [mM]	Comment	Tit. [µl]	A	B
			~200 µM O2		MiR			
AmR			10 µM	10		2		
HRP			1 U/ml		500 U/ml stock	4		
SOD	OBL		5 U/ml		5000 U/ml stock	2		
		1HP0	0		<i>Mark on Amp signal</i>			
1HP.1		1HP.1	0.1 µM	0.04	<i>Mark on Amp signal</i>	5		
1HP.2		1HP.2	0.2 µM	0.04	<i>Mark on Amp signal</i>	5		
mt	ROX							
		2HP0	0		<i>Mark on Amp signal</i>			
2HP.1		2HP.1	0.1 µM	0.04	<i>Mark on Amp signal</i>	5		
1P			5 mM	2000		5		
1M	1PM		2 mM	400		10		
2D	2D		2.5 mM	500		10		
3U	3U		Δ0.5 µM	1	CCCP	Δ1		
4G	4G		10 mM	2000		10		
		3HP0	0		<i>Mark on Amp signal</i>			
3HP.1		3HP.1	0.1 µM	0.04	<i>Mark on Amp signal</i>	5		
5S	5S		50 mM	1000		100		
6Oct	6Oct		0.5 mM	100		10		
7Rot	7Rot		0.5 µM	1		1		
		4HP0	0		<i>Mark on Amp signal</i>			
4HP.1		4HP.1	0.1 µM	0.04	<i>Mark on Amp signal</i>	5		
8Gp	8Gp		10 mM	1000		20		
9Ama	9Ama		2.5 µM	5		1		
		5HP0	0		<i>Mark on Amp signal</i>			
5HP.1		5HP.1	0.1 µM	0.04	<i>Mark on Amp signal</i>	5		