


FNSGp(PGM) 1R;1Dig;1D;1M.1;2Oct;2c;2M2;3P;4G;5S;6Gp;7U;8Rot;9Ama;10Tm;11Azd



O2k and DatLab file: P___ (A / B)		2017-		Operator:				
Sample type:		Cohort:		Sample code:				
Sample.Subsample number:		Unit:		Concentration:				
Medium: MiR06Cr								
Event	Mark name	State	Final conc. 2 ml O2k	Stock [mM]	Comment	Tit. [μl]	A	B
			~200 μM O2		MiR			
1ce	1R	R						
1Dig	1ROX	ROX		8.1				
1D	1D	ROX	1 / 2.5 mM	500		4 / 10		
1M.1	1M.1		0.1 mM	50		4		
2Oct	2Oct	F _P	0.5 mM	100		10		
2c	2c	F(N) _{PC}	10 μM	4		5		
2M2	2M2	F(N) _P	2 mM	400		9.5		
3P	3P	FN(PM) _P	5 mM	2000		5		
4G	4G	FN(PGM) _P	10 mM	2000		10		
5S	5S	FNS(PGM) _P	50 mM	1000		100		
6Gp	6Gp	FNSGp(PGM) _P	10 mM	1000		20		
7U*	7U	FNSGp(PGM) _E	Δ0.5 μM	1	CCCP	Δ1		
8Rot	8Rot	SGp _E	0.5 μM	1		1		
9Ama	9Ama	ROX	2.5 μM	5		1		
10As			2 mM	800		5		
10Tm			0.5 mM	200	Immediately after As	5		
10O					~20 min open			
10C	10Tm	Tm _E			Close and wait 5 min			
11Azd	11Azd	ROX	≥100 mM	4000	~5 min	100		